

Dosing of ICG for Xiralite			
1 vial 25 mg ICG, add 5 ml water for injection (equals 5 mg / ml)			
Bodyweight (lbs)	Bodyweight (kg)	mg ICG at 0.1 mg/kg BW	Amount of ICG (ml)
44	20	2	0.4
55	25	2.5	0.5
66	30	3	0.6
77	35	3.5	0.7
88	40	4	0.8
99	45	4.5	0.9
110	50	5	1.0
121	55	5.5	1.1
132	60	6	1.2
143	65	6.5	1.3
154	70	7	1.4
165	75	7.5	1.5
176	80	8	1.6
187	85	8.5	1.7
198	90	9	1.8
209	95	9.5	1.9
220	100	10	2.0
231	105	10.5	2.1
243	110	11	2.2
254	115	11.5	2.3
265	120	12	2.4
276	125	12.5	2.5
287	130	13	2.6
298	135	13.5	2.7
309	140	14	2.8
320	145	14.5	2.9
331	150	15	3.0